


# LC8010 Isocratic Liquid Chromatography System

## LC8010 | Chromatography Systems

Isocratic HPLC system for stable routine methods where a fixed mobile-phase composition is required.



### Key Features

- Isocratic liquid chromatography workflow
- Suitable for routine laboratory methods and teaching laboratories
- Image set includes product, detail, and laboratory-background photos
- Can be quoted as a compact HPLC configuration

### Typical Applications

- Routine QC methods
- Teaching laboratories
- Food and environmental analysis
- Method validation with isocratic conditions

### Selection Reference

|                        |  |
|------------------------|--|
| <b>Instrument Type</b> | Isocratic HPLC   |
| <b>Separation Mode</b> | Fixed mobile-phase composition                                   |
| <b>Typical Output</b>  | Chromatographic separation and quantitative analysis             |
| <b>Configuration</b>   | Pump, injector, column, detector, and workstation by requirement |
| <b>Voltage</b>         | Confirm voltage and plug standard before quotation               |
| <b>Document Status</b> | English PDF version for international product selection          |

Final configuration, accessories, voltage, plug standard, documentation scope, and delivery details should be confirmed before quotation.

## Product Views

